

IsoAir® 310P

Aerosol Particle Sensor



The IsoAir 310P Aerosol Particle Sensor makes aerosol monitoring easy and cost-effective. It has unparalleled performance in a rugged, chemical-resistant, and easy-to-disinfect 316L stainless steel box. The clean, reliable internal pump lowers installation costs and the need for costly vacuum system connections.

FEATURES

- 4 standard channel sizes with 2-channel Pharmaceutical mode
- Flow rate of 1.0 CFM
- Passivated 316L stainless steel enclosure
- Pump with HEPA filtered exhaust
- Automatic shut-off for capped probe
- Compliance with ISO 21501-4 and CE standards
- Two local LEDs for status and activity
- Optional 4-20 mA output
- Ethernet, MODBUS, or OPC communications
- Optional 4-20 mA input and solid state contact outputs for remote status indication



BENEFITS

Ensure Compliance

- Real-time particle monitoring gives quick notification of particle anomalies
- Queues 3,000 samples in the event of computer system failure
- Detects 0.5 and 5.0 μm for GMP and FDA compliance
- Ethernet and 4-20 mA outputs meet 21 CFR Part 11 regulations for data transmission security
- Identifies breakdown in isolator protection

Increase Productivity

- Interfaces with optional Pharmaceutical Net or Facility Net software for comprehensive management of the cleanroom
- Automatic paging allows rapid response to alarms
- Optional 4-20 mA inputs allow for integration of up to 3 environmental input sensors
- Easy to remove sensor simplifies calibration and service
- Validation documentation available

Cost-Effective

- Complete monitor is easy to install
- No central vacuum system needed
- Rugged, chemically-resistant, easy to disinfect enclosure
- Passivated stainless steel surface has well defined microbial control characteristics
- Diode laser reduces maintenance requirements

APPLICATIONS

- Pharmaceutical manufacturing area
- Isolator monitoring
- Cleanroom monitoring
- Aerospace and industrial applications
- Trend analysis
- Statistical process control

Without measurement there is no control.

IsoAir® 310P

Aerosol Particle Sensor

specifications

Channels	0.3, 0.5, 1.0 and 5.0 µm channels or can be set for two-channel mode: 0.5 and 5.0 µm
Flow rate	1.0 ft³/min (28.3 L/min)
Sample probe/tubing:	3/8" ID; up to two meters of 3/8 inch tubing can be connected to the inlet
Vacuum source	Internal reciprocal blower is cool running long-lived with no contaminating carbon vanes. Blower output is HEPA filtered with exhaust out bottom of enclosure.
Blower shut off	Automatic if blower fails or if system detects sample probe capped while pump is on
Calibration	Calibration materials used are traceable to the US National Institute of Standards and Technology. Unit is calibrated to ISO 21501-4.
Zero count level	≤ 7.07 counts/m³
Maximum concentration @ 10% coincidence loss	1,374,269/ft³
Exterior surface	Passivated stainless steel (Grade 316L)
Communications	Ethernet (data, instructions); RS-232 (set-up and diagnostics only) Modbus TCP; OPC communications Optional 4-20 mA output (2 data channels plus 1 status channel) Optional 5 solid-state outputs and 3 4-20 mA inputs
Controlling software	Pharmaceutical Net, Facility Net (must purchase separately)
Data security	Unalterable data; 21 CFR Part 11 compatible transfer via Pharmaceutical Net. Backup storage of 3,000 samples.
Status indicators	Ethernet Programmable Programmable: 3-color LED Activity: 1-color LED 4-20 mA Interface Laser and flow status: 3-color LED Activity: 1-color LED
Dimensions (w, h, d)	10 x 12 x 6 in (25 x 30 x 15 cm)
Weight	20 lb (9 kg)
Power	100 – 240 VAC, 1 A, 50/60 Hz universal power supply
Environment	Temperature: 39 – 95 °F (4 – 35 °C) Humidity: 5% – 95% Non-condensing Altitude: 0 – 6,562 ft (0 – 2,000 m)
Standard	Sample probe (stainless steel), probe cap and chain; spare fuse(s); tubing
Manual	IsoAir 310P Operations Manual included
Zero count filter	Included

IsoAir® is a registered trademark of Particle Measuring Systems, Inc.
Particle Measuring Systems, Inc. reserves the right to change specifications without notice.



Particle Measuring Systems Headquarters
5475 Airport Blvd., Boulder, CO 80301, USA
Tel: +1 303 443-7100 +1 800 238-1801

FAX: +1 303 449-6870

Instrument Service & Support: +1 800-557-6363

Customer Response Center: +1 877-475-3317

PMS EMEA
Tel: +44 1684 581 000
Email: PMSEMEA@pmeasuring.com

Particle Measuring Systems France
Tel: +33 160 10 32 96
Email: PMSFrance@pmeasuring.com

Particle Measuring Systems Germany
Tel: +49 2222 9299 522
Email: PMSGermany@pmeasuring.com

Particle Measuring Systems Italia
Tel: +39 06 9053 0130
Email: PMSIRL@pmeasuring.com

Particle Measuring Systems Nordic
Tel: +45 707 028 55
Email: PMSNordic@pmeasuring.com

Particle Measuring Systems China
Tel: +86 21 6113 3600
Email: pmschina@pmeasuring.com

Particle Measuring Systems Japan
Tel: +81 3 5298 8175
Email: PMSJapan@pmeasuring.com

Particle Measuring Systems Singapore
Tel: +65 6496 0330
Email: PMSSingapore@pmeasuring.com

AUTHORIZED REPRESENTATIVE

Particle Measuring Systems Brazil
Tel: +55 11 5188 8166
Email: PMSBrazil@pmeasuring.com

Particle Measuring Systems Mexico
Tel: +52 55 2271 5106
Email: PMSMexico@pmeasuring.com

Particle Measuring Systems Puerto Rico
Tel: +1 787 718 9096
Email: PMSPuertoRico@pmeasuring.com